



Microdacyn® Wound Care Solution

Microdacyn® Hydrogel

Super-oxidized hypochlorous acid solutions for advanced tissue repair



Antimicrobial | Biofilm | Healing | How to Use | Clinical Evidence

Microdacyn is a class IIb Medical Device for use in the debridement and moistening of acute and chronic wounds, ulcers, cuts, abrasions and burns including those located in any human cavity such as the oral, nasal or ear. Sodium hypochlorite and hypochlorous acid are ancillary substances may have a local antimicrobial effect. Through reducing the microbial load and assisting in creating a moist environment, it enables the body to perform its own healing process. Microdacyn® Wound Care can be broadly applied within a comprehensive wound treatment. Do not use if sensitive to hypochlorous acid or sodium hypochlorite.

TeArai
BioFarma


Antimicrobial | Biofilm

Super-oxidized solutions (SOS) and hydrogels utilize physiological concentrations of reactive oxygen species such as hypochlorous acid to induce osmolysis in a totally non-antibiotic & chemical free process

Hypochlorous acid is naturally present in the human body & is produced as a natural response to infection by the white blood cells in our immune system such as neutrophils through a mechanism known as an oxidative burst.

Unlike traditional antiseptics & wound washes Microdacyn® only causes damage to single cell microbes such as bacteria and is not cytotoxic to human cells allowing the delicate cells involved in healing to progress unhindered¹⁵

Recent consensus²⁸ has concluded that **rather than inhibit wound healing, Microdacyn® provides improvement to wound healing** and is indicated for use without restriction in delicate tissue due to its low cytotoxicity.



Proven effective against: Bacteria (incl MRSA), Fungi, Viruses, Spores & Uniquely Penetrates Biofilms²

“Use antiseptics at the lowest effective concentration to minimize harm to skin and tissue cells involved in wound healing”²¹

International Wound Infection Institutes (IWII) Clinical Practice Guidelines 2016


Treatment	Betadine®	Prontosan®	Microdacyn®
Generic Name	Antiseptic	Surfactant-based Antiseptic	Super-Oxidised Solution
Concentration	10%	0.1%	0.004%
Mechanism	Chemical	AMP	Osmolysis

The enormous difference in germicidal potency is due to the fact that pure HOCl as an uncharged species can penetrate microbial cell walls, whereas charged ionic species cannot¹


Super-Oxidised Solutions Do Not Promote Bacterial Resistance²¹

Biofilm

Super-oxidised solutions *“rapidly penetrate biofilms killing microbes from within”²¹*



t=0



t=10

An example of this rapid penetration is shown here:

Fluorescence microscopy of a mature, 6 day, *Pseudomonas aeruginosa* biofilm before and after 10 minutes of repeated exposure. It is clear to see uniform penetration and eradication of living microbes

Green | Living Microbes Red | Dead Microbes

Inflammation | Healing | Clinical Evidence



The low concentration Hypochlorous acid found in Microdacyn® is the only antimicrobial component with clinical research demonstrating additional wound healing effects.


The use of Microdacyn has been demonstrated to:

- **Reduce Inflammation, Itch & Irritation**²⁰
- **Increase available oxygen (TcPO2 - mmHg)**^{3,5}

Here we observe the average tissue oxygenation effects when applying Microdacyn® to 393 non-insulin dependent diabetic (NIDDM) patients velous leg ulcers.

The TcPO2 levels are significantly increased to quickly improve tissue conditions within the wound and are sustained above baseline for up to 72 hours⁵

Effect of comorbidities on microcirculatory response of periwound tissue to Microdacyn® exposure.



These additional effects may further assist in Microdacyn® accelerating wound healing & being excellently tolerated

Clinical Evidence

Clinical conclusions of 10+ years of Microdacyn use internationally



Improves Healing

- **Improved Healing**
 - Improve wound healing^{5-9,11,13-15}
 - Reduce wound healing time^{5-9,11,13-15}
 - Improve wound bed oxygen supply (TcPO2)^{3,5}
 - Improve wound bed granulation^{3,11,13}
 - To be non-cytotoxic^{1-2,11,14,19,21}



Prevents Infection

- **Prevents Infection**
 - Improve infection control^{2,4,7,14,16}
 - Reduce the requirement for antibiotic use to manage infection¹⁰
 - Reduce post-surgical infection rate^{2,7,12,15}



Reduces Inflammation

- **Patient Well Being**
 - Reduce pain associated with wound cleansing¹⁷
 - Reduced wound associated malodour¹¹
 - Reduce patient length of stay^{8,16-17}
- **Additionally**
 - To be at least as effective as certain oral antibiotics for treatment of mild diabetic foot infections¹⁰
 - To be a wider spectrum antimicrobial compared to commonly used topical-antimicrobials^{12,17-18}

How to Use | Product Details


Microdacyn® should be applied liberally at every dressing change or if the wound can be regularly accessed, up to 3 times daily

1 Saturate



Irrigate the wound by applying Microdacyn® Wound Care Solution directly onto/into the wound to the point of saturation

2 Soak



Place Microdacyn® saturated gauze onto the wound for up to 10 minutes helping to hydrate and aid slough removal

Microdacyn® soaks are to assist manual cleaning/debridement not to replace it

3 Flush Clean



Once the wound is clean use Microdacyn® Wound Care Solution is an irrigant to flush any remaining debris clear prior to dressing

Practice Tip

Another strip of gauze or similar that is doused in Microdacyn® Wound Care Solution is excellent for cleaning both in & around the wound

4

Cover & Dress

Leave Microdacyn® in the wound, do not rinse out

3b



Higher Risk Patients

In the case of the following, repeat step 2;

- Significant Debridement
- Chronic Wounds
- Non-Healing Wounds
- Suspected Biofilm

3c



Applying Hydrogel

Leave Microdacyn® Wound Care Solution in the wound and apply Microdacyn® Hydrogel 2-5mm thick as required
This may be left for up to 3 days
Microdacyn® Hydrogel is well suited to painful, inflamed or malodorous wounds

Microdacyn® Wound Care

Super-Oxidised Water, Sodium Chloride (0.022%), Hypochlorous Acid (0.004%), Sodium Hypochlorite (0.004%), Ancillary Substances including Oxygen & Ozone

Microdacyn® Hydrogel

Super-Oxidised Water, Sodium Chloride (0.066%), Hypochlorous Acid (0.004%), Sodium Hypochlorite (0.002%), Ancillary Substances including Oxygen & Ozone
Additionally Sodium Magnesium Fluorosilicate Gelling Agent (3.0%)

Product Description	Internal Code	
Microdacyn® Wound Care Solution 120ml	MDWC120	Box 24
Microdacyn® Wound Care Solution 250ml	MDWC250	Box 12
Microdacyn® Surgical Irrigation Solution 990ml	MDSIWT990	Box 6
Microdacyn® Hydrogel 60g	MDHG60	Box 24
Microdacyn® Hydrogel 120g	MDHG120	Box 24

**“Stop anointing wounds
& start cleansing wounds”**

International Wound Infection Institutes
Clinical Practice Guidelines

In the case of severe or multi-resistant infection, stubborn biofilm or lack of progress.
Microdacyn® may be applied copiously to achieve desired results

Microdacyn® is compatible with:

- All Wound Dressings
- Leaving Soaked Gauze in the Wound
- Pulse Lavage / Wound Irrigation
- Under Occlusion
- NPWTi & Ultrasonic Debridement

Microdacyn® can be used on:

- Skin
- Mucosa
- Exposed Ligament, Joints, Bones & Tendons
- Sensitive Areas Such as Eyes
- Children
- ANY Human Cavity

Product Information

- Ready-to-use, pH neutral solution
- 24 month unopened shelf life
- Discard Within 60 Days of Opening
- Does not promote bacterial resistance²¹
- Available from your local pharmacy

www.microdacyn.com.au | 1800 Te Arai (83 2724) | enquiries@tearaibiofarma.com



Microdacyn®
.com.au